

## Foaming Flux Remover in 30ml (1.0oz) Squeeze Bottle w/tip

### Product Highlights

Foaming Flux Remover

Can be used with Leaded and Lead-Free applications

**RoHS 3 and REACH compliant**

### Specifications

Packaging:

30ml Squeeze Bottle

Shelf Life:

Refrigerated >24 months, Unrefrigerated >24 months



### Storage and Handling

Store refrigerated or at room temperature 3-25°C (37-77°F). The flux remover can be allowed to freeze. Freezing will not degrade this product. Allow 4 hours for flux remover to reach an operating temperature of 20-25°C (68-77°F) before use.

### Transportation

This product requires ground shipping. Shipping below 0°C (32°F) or above 25°C (77°F) for normal transit times by ground will not impact this product's stated shelf life.

### Conforms to the following Industry Standards:

RoHS 3 Directive (EU) 2015/863:

Yes

## Flux Remover Selection Guide

Description	Packaging	Orderable Part Number
Degreaser / Contact Cleaner (Non-Flammable)	Pen, 10ml	CQNRM
	Squeeze Bottle, 0.5oz	CQNRM-0.5
	Squeeze Bottle, 1.0oz	CQNRM-1.0
Flux Remover (Non-Flammable)	Pen, 10ml	CQNTRM
	Squeeze Bottle, 0.5oz	CQNTRM-0.5
	Squeeze Bottle, 1.0oz	CQNTRM-1.0
Flux Remover (Fast Evaporating, Non-Flammable)	Pen, 10ml	CQNIRM
	Squeeze Bottle, 0.5oz	CQNIRM-0.5
	Squeeze Bottle, 1.0oz	CQNIRM-1.0
Flux Remover (Foam-Free, Flammable)	Pen, 10ml	CQIRM
	Squeeze Bottle, 0.5oz	CQIRM-0.5
	Squeeze Bottle, 1.0oz	CQIRM-1.0
Flux Remover (Foaming, Flammable)	Pen, 10ml	CQ6RM
	Squeeze Bottle, 0.5oz	CQ6RM-0.5
	Squeeze Bottle, 1.0oz	CQ6RM-1.0